ROLLY
890.01

CE
Regolamento (UE)
2017/745

IT

WWW.KONG.IT

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INDICE

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Download the digital version

TABLE OF CONTENTS

1 - SYMBOLS AND SUPPORT 4
1.1 Symbols 4
1.2 Support 4

2 - GENERAL INFORMATION 5

3 - TECHNICAL FEATURES 6
3.1 Terminology and materials of the parts 6
3.2 Technical data 8
3.3 Optional components and spare parts 9

4 - SPECIFIC INFORMATION 10
4.1 Intended use 10
4.2 Preparation 10
4.3 Setup 11
4.4 Blocking the patient 16
4.5 Transporting the patient 18

5 - MAINTENANCE AND REPAIRS 20
5.1 General 20
5.2 Maintenance 20
5.3 Repair 20

6 - STORAGE 21

7 - CHECKS, INSPECTIONS AND SERVICING 22
7.1 Checks 22
7.2 Inspections 22
7.3 Servicing 22

8 - PRODUCT LIFE AND GUARANTEE 23
8.1 Product life 23
8.2 Disposal 23
8.3 Guarantee 23
8.4 Law obligations 23

9 - LABELLING AND SYMBOLS 24
9.1 Medical device label 24
9.2 Symbols 24

10 - DOCUMENTS 25
10.1 Inspections and servicing register 25
10.2 Maintenance and repair register 26
10.3 EU Declaration of conformity 27

11 - LIST OF APPLIED STANDARDS AND REFERENCE 28
11.1 Applied standards 28
11.2 Standards used as reference 29

CHAPTER 1

SYMBOLS AND SUPPORT

1.1 SYMBOLS

To make reading this manual comfortable and clear, the following are the symbols used to manage important warnings, for correct and safe use of the device.

REQUISITE FOR CORRECT USE
It identifies the presence of information required to use the device correctly.

INFORMATIVE REQUISITE
It identifies the presence of useful and general information, to guide the user to consciously use the device and/or carry out the actions.

It identifies that the Medical Device is manufactured, designed and produced in compliance with the provisions of the General Safety and Performance Requirements of Regulation (EU) 2017/745 (Class I medical device, in compliance with classification rule 1 as indicated in the attachment VIII).

1.2 SUPPORT

For information, contact Kong Customer Service by:
- telephone +39 0341 630506,
- fax +39 0341 641550,
- email: safetycare@kong.it,

or write to KONG S.p.A. – Via XXV Aprile, 4 – 23804 Monte Marenzo LC - ITALY.
To facilitate assistance operations, always communicate or indicate the serial number (SN) indicated on the label applied to the Medical Device.

CHAPTER 2

GENERAL INFORMATION

The information provided by the manufacturer (hereinafter referred to as information) must be read and clearly understood by the user prior to using the Medical Device (MD hereafter). The information regards the description of the features, performance, assembly, disassembly, maintenance, preservation, disinfection, etc. of the device. Even though the information offers tips on use, this information shall not be considered as a user manual under actual conditions of use.

2.1 WARNINGS AND LIMITATIONS OF USE:

- this device shall be strictly used by people who are physically fit, trained (instructed and taught) to use the device and with specific experience regarding moving the patient or, during training activities, by people under direct supervision of the trainers/ supervisors who guarantee the safety thereof,
- do not use the device before fully reading and understanding this user manual,
- strictly follow the manufacturer's information, improper use of the device is hazardous,
- modifying and/or repairing the device is strictly forbidden,
- all checks described in chapter 7 must be carried out prior to and after using the device. In case of any doubt on the efficiency of the device, the user must replace it immediately,
- non-compliant use, deformations, falls, wear, chemical contamination, exposing textile/ plastic components/devices to temperature below -30°C or above +50°C and metal components/devices to temperatures exceeding 100°C are some examples of causes that can reduce, limit and end the life of the device,
- Before use, ensure that the weight does not exceed the load capacity tested for this MD,
- in order to reduce risks of exposure to / transmission of infectious diseases, clean and disinfect the device as indicated in chapter 5,
- improper use of the patient immobilisation systems can jeopardise the safety of the patient,
- always check the compatibility of the devices used alongside the device by consulting the relative manufacturer's information,
- use of spare parts or optional components different from the ones defined for this MD,
- do not expose the device to sources of heat and at contact with chemical substances. Reduce exposure to direct sunlight as much as possible. At low temperatures and in the presence of humidity, formation of ice could reduce flexibility and increase cutting and abrasion-related risks on textile and synthetic devices,
- report to the manufacturrer and the competent authority of the member state where the user and / or patient is established any serious incident occurring in relation to the MD.

All our devices are tested/ checked piece by piece in compliance with the procedures laid down by the Quality System certified in accordance with the UNI EN ISO 9001 standard. Laboratory tests, testing, information and standards do not always reproduce the practical result. Thus, the results obtained under the actual conditions of use of the device may differ, even considerably at times. The best information lies in the continuous practical use under the supervision of skilled/expert/qualified people.

⚠ **Warning: not suitable for use in an ATEX environment (Directive 94/9/CE)**

CHAPTER 3

TECHNICAL FEATURES

3.1 TERMINOLOGY AND MATERIALS OF THE PARTS

A - Support surface in high-density polyethylene
B1 - Slings in polyester
B2 - Slings in polyester
C - Slings in polyester with eyelets
D - Slings in polyester
E - Sling in polyester
F - Handles in polyester
G - Hanging kit in high-tenacity polyester (G1 short part - G2 long part)
H - Static rope
I - Buckle with hook in aluminium alloy
L - Buckle with hook in aluminium alloy
M - Slings with eyelets
N - Sling with ratchet
O - Transport bag.


Fig.1

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
3.2 TECHNICAL DATA

3.2.1 Weight and dimensions

Length:	245 cm
Width:	92 cm
Rolled in the bag:	Ø 30 x 110 cm
Overall weight of the parts (fig. 1):	9 kg

3.2.2 Tested load capacities

Mode of use	**Tested load capacities**
Transport with the handles (F) (fig. 16)	250 kg
Lifting and lowering with the hanging kit (G) (fig. 17)	1500 kg
Vertical lowering and recovery (fig. 18)	1200 kg

3.2.3 Buoyancy force

Buoyancy force of the kit 880020300KK FLOK-IT (optional component)	400 N (~40 kg)

Safety factors shall be defined by the user.

3.3 OPTIONAL COMPONENTS AND SPARE PARTS

The methods of connection and use of the optional components are defined in the relative user instructions.

3.3.1 Optional components

880020300KK FLOK IT
Floating kit for stretchers, made of three cylinders of different lenghts (2 side pieces of 1 m each and 1 piece head side of 0,4 m) , each of one equipped with two fastening straps with quick release plastick buckles.


Fig.3

3.3.2 Spare Parts

B1 - Slings in polyester
B2 - Slings in polyester
C - Slings in polyester with eyelets
D - Slings in polyester
E - Sling in polyester
F - Handles in polyester
G - Hanging kit in polyester high tenacity (G1 short part - G2 long part)
H - Static rope
N - Sling with ratchet
O - Transport bag.


Fig.4

CHAPTER 4

SPECIFIC INFORMATION

4.1 INTENDED USE

The Mdical Device (MD) **“ROLLY”** is a rescue sheet that is particularly suited to recovering and transporting a patient - even immobilised on **"X-TRIM"** spinal boards or on the **"VACUUM”** mattresses, in a quick manner and in particularly difficult situations.
Decisions concerning moving and immobilising the patient, as well as the duration, method to be used and combination with other devices, shall be taken and executed by expert and trained personnel only.
The supplied hanging kit (G) makes it possible to use the “ROLLY” sheet with a winch.

4.2 PREPARATION

a) Remove the "ROLLY" rescue sheet from the transport bag (O) and remove the ratchet sling (N)
b) Distend the support surface (A) and flatten it, rolling the ends over in the opposite direction and folding the central section (fig. 5)


Fig. 5

4.3.1 Setting up for blocking the patient

Prepare the following setup to block the patient:

a) Insert the sling (E) in the A4 holes in the support surface (fig. 6).
b) Insert the slings (B1 and B2) in the A1 eyelets in the support surface (fig. 7).
c) Insert the slings with eyelets (C) in the A2 eyelets in the support surface (fig. 8).
d) Insert the slings (D) in the A3 eyelets in the support surface (fig. 9).


Fig.6


Fig.7


Fig.8



Fig.9

4.3.2 Setting up for moving the patient

In relation to the lifting and transport procedures, which the rescue workers considers suitable, set the stretcher up with at least one of the following parts:
a) manual transport: insert only the side handles (F) in the eyelets A5 on the support surface (fig. 10) - Please note: we suggest inserting the shortest couple of handles in the place near the torso of the patient, in order to keep the head higher than the feet.


Fig.10

b) lifting/lowering and/or transportation using the lifting kit (G): insert the short hanging kit (G1) in the eyelets A6 and the long hanging kit (G2) in the eyelets A7 on the support surface (fig. 11).


Fig.11

c) lifting/lowering or vertical lowering: insert the rope (H) in the holes A8 in the support surface (fig. 12).


Fig.12

4.4 BLOCKING THE PATIENT

After having positioned the patient on the support surface (A):

a) Insert the buckles with hook (I) in the corresponding slings with eyelet (C) and tension the slings (B) to lift the side parts of the support surface (A) enough to contain the patient (fig. 13).

b) Insert the buckles with hook (L) in the corresponding slings with eyelet (M) and tension the slings (D) to lift the end part of the support surface (A) enough to contain the patient's feet (fig. 14).


Fig.13


Fig.14

c) Tension the sling (E) to lift the end part of the support surface (A) enough to contain the head of the patient and immobilise it (fig. 15).
Avoid prolonged contact of the patient's skin with the fabric of the device.

This procedure must also be carried out with patients immobilised on the spinal board "X-TRIM" or on the vacuum-operated mattress "VACUUM".


Fig.15

4.5 TRANSPORTING THE PATIENT

The MD **“ROLLY”** is suitable for use in the following modes of use:

- **Manual transportation (fig. 16)**
 with the handles (F)

- **Lifting and lowering (fig. 17)**
 with the hanging kit (G)

- **Vertical descent and recovery (fig. 18)**
 with the rope (H)

- **Water transport (fig. 19)**
 with the kit 880020300KK FLOK-IT (P)

Before any recovery operation, make sure that the weight does not exceed the tested load capacity!


Fig.16


Fig.17


Fig.18


Fig.19

CHAPTER 5

MAINTENANCE AND REPAIR

5.1 GENERAL INFORMATION

This MD is made of material which is highly resistant to wear and external agents. Despite this, the use conditions require maintenance operations and even repair in particular conditions.

Maintenance and repair operations must be recorded in a special register, an example of which is indicated in chapter 10.

5.2 MAINTENANCE

The maintenance operations that must be carried out by the user are cleaning and disinfection (when necessary).

Cleaning: wash the device by hand in lukewarm water (30-32°C) with neutral detergent (pH 5.5-8.5) and allow it to air dry, avoiding direct exposure to the sun and heat sources,

Disinfection: immerse the MD (except for devices in/with parts in Dyneema® and aramid e.g. Kevlar®) for at least 30 minutes in warm water (58-60°C) and allow it to air dry, avoiding direct exposure to the sun and heat sources. Do not bleach, do not wash in a washing machine, do not spin-dry, do not tumble dry, do not iron and avoid sterilizing textile devices in an autoclave.

5.3 REPAIR

Repair shall be strictly carried out by the manufacturer.
The user is allowed to replace the parts mentioned in paragraph 3.3.2 only, with new and genuine parts.

CHAPTER 6

6 STORAGE

After cleaning, disinfecting (if necessary) and drying, store the MD and its accessories in a dry (relative humidity 40-90%), cool (temperature 5-30°C), dark, well-ventilated place, away from sharp edges, chemicals and other harmful conditions.

Do not store this device wet!

CHAPTER 7

CHECKS, INSPECTIONS AND REVISIONS

7.1 CHECKS

Before and after every use, check the device and ensure that:
- it is appropriate for the use it was intended for,
- it has not been subjected to mechanical deformations and it neither has cracks nor shows signs of wear,
- the textile parts do not have cuts, burns, chemical products residues, excessive hair, wear, in particular check areas at contact with the metal components (buckles, rings, etc.);
- the markings, including labels, are legible.

7.2 INSPECTIONS

The device must be inspected at least once a year, starting from the date of first use, by expert personnel trained and approved by KONG S.p.A. The date of the first use and the outcome of the inspection must be recorded in the Inspections and Servicing Register, regarding which an example is indicated in chapter 10.

7.3 REVISIONS

The device must be serviced by KONG S.p.A., or by expert personnel trained and approved by KONG S.p.A. if:
- any malfunctions are detected,
- the result of the pre and post use check or of the inspections is negative.

The outcome of the servicing must be recorded in the Inspections and Servicing Register, regarding which an example is indicated in chapter 10.
The revised device is guaranteed for one year from the revision date.

CHAPTER 8

DURATION AND WARRANTY

8.1 DURATION

Carefully read the “Warnings and limitations of use” paragraph in chapter 2.
The duration of this device depends on the outcome of the at least one-yearly inspection (paragraph 7.2) and possible servicing.
The hanging kit, even if not used, must be replaced after 10 years from the first use and in any case no later than 12 years from the date of manufacture.

Important: any devices that have not passed the pre-use, post-use and periodic inspection must be eliminated and made unusable.

8.2 DISPOSAL

Follow the rules in force in the country of use and the waste disposal procedures of the reference hospital facility for proper disposal.

8.3 WARRANTY

The manufacturer guarantees the compliance of the device with the regulations in force at the time of production. Warranty covering defects is limited to defects regarding raw materials and manufacturing defects. The warranty does not cover normal wear, oxidation, damage caused by non-compliant use and/or in competition, improper maintenance, transportation, preservation or storage. Warranty will be immediately deemed null and void should the device be modified or tampered with. Validity corresponds to the legal validity of the country where the device was sold, as from the date of sale, by the manufacturer. No claims shall be raised against the manufacturer once the aforementioned period has expired. Any repair or replacement request under warranty must be accompanied by a proof of purchase. Should the defect be recognised, the manufacturer undertakes to repair or, replace or refund the device at his discretion. The manufacturer shall not be held liable beyond the invoice price of the device in any case whatsoever.


8.4 LAW OBLIGATIONS

Professional and free time activities are often regulated by special national laws which can impose limits and/or obligations concerning the use of these devices. The user has the obligation to know and apply such laws which could provide for limits different from those indicated in this information.

CHAPTER 9

LABELLING AND SYMBOLS

9.1 MEDICAL DEVICE LABEL

Any changes to the positioning of the symbols does not change the content thereof.

9.2 SYMBOLS

Identification and traceability

MD	Medical Device
REF	Product identification code
UDI	Unique Device Identifier
	Manufacurer's identification
SN	Unique serial number
ITALY YYYY / MM	Country, year and month of production.

Data matrix UDI

Information and safety

	Consult the user manual
CE	Compliant with Regulation (EU) 2017/745

CHAPTER 10

DOCUMENTS

10.1 INSPECTIONS AND SERVICING REGISTER

RESCUE SHEET - Class I MD			
REF.	890.01	**MOD.**	ROLLY
UDI-DI	08023577051727	**SN**	
Date of first use			

INSPECTION AND SERVICING				
DATE	**I/S**	**Description**	**Outcome**	**Supervisor**

10.2 MAINTENANCE AND REPAIR REGISTER

RESCUE SHEET - Class I MD			
REF.	890.01	**MOD.**	ROLLY
UDI-DI	08023577051727	**SN**	
Date of first use			

MAINTENANCE AND REPAIRS

DATE	**M/R**	**Description**	**Outcome**	**Supervisor**

10.3 EU DECLARATION OF CONFORMITY

The EU Declaration of Conformity is available for download at www.kong.it

KONG s.p.a.
via XXV Aprile, 4 | 23804 Monte Marenzo (LC) - ITALY | tel +39 0341.630506 | info@kong.it

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www.kong.it

CHAPTER 11

LIST OF APPLIED AND REFERENCE STANDARDS

11.1 APPLIED STANDARDS

- EN ISO 14971: 2019
 Medical devices - Application of risk management to medical devices.
- EN 62366-1: 2015 + A1: 2020
 Medical devices - Application of usability engineering to medical devices.
- EN ISO 15223-1: 2021
 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements
- EN ISO 20417: 2021
 Medical devices - Information to be supplied by the manufacturer
- EN ISO 10993-1: 2020
 Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
- EN ISO 10993-5: 2009
 Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
- EN ISO 10993-23: 2021
 Biological evaluation of medical devices - Part 23: Tests for irritation
- EN 597-1: 2015
 Furniture - Assessment of the ignitability of mattresses and upholstered bed bases - Part 1: Ignition source smouldering cigarette
- EN ISO 846:2019
 Plastics - Evaluation of the action of microorganisms
- EN 1865-1: 2010 + A1: 2015
 Patient handling equipment used in road ambulances - Part 1: General stretcher systems and patient handling equipment

11.2 STANDARDS USED AS REFERENCE

- EN ISO 13485: 2016
 Medical devices - Quality management systems - Requirements for regulatory purposes
- EN 13718-1: 2014 + A1: 2020
 Medical vehicles and their equipment - Air ambulances - Part 1: Requirements for medical devices used in air ambulances.